



















Components description:


Name	Description
<p data-bbox="304 264 643 293">48-tube magnetic homogenizer</p> 	<p data-bbox="735 271 1473 439">48-tube magnetic homogenizer is installed on the worktable of the device in versions DTstream9 L1, DTstream9 L4, DTstream12 L1, DTstream12 L4, DTstream15 L1, DTstream15 L4 and is designed for the homogenization of the contents of test tubes by rotating a magnetic rod placed in the tube, under the influence of the magnetic field of the rotor located in the body of the homogenizer.</p> <p data-bbox="735 439 1473 521">The rotor is rotated by a micro electric motor mounted in the body frame. To power the drive micro-motor and control the operation, magnetic homogenizer is connected to the CON socket of the device.</p> <p data-bbox="735 521 1473 577">Dimensions of the 48-tube magnetic homogenizer: 133±0,2 x 220±0,2 x 78±0,2 mm.</p> <p data-bbox="735 577 1473 607">Weight of the 48-tube magnetic homogenizer: 1180±2 g.</p> <p data-bbox="735 607 1473 663">Structurally the magnetic homogenizer consists of a base (body), with a magnetic rotor with a drive inside of it.</p> <p data-bbox="735 663 1473 714">48-tube magnetic homogenizer body and rack are made of Trovidur EN grade polyvinyl chloride.</p> <p data-bbox="735 714 1473 797">48-tube magnetic homogenizer is installed on the worktable of the device and fixed on the worktable by aligning the pins located on the worktable of the device and holes in the base of a homogenizer.</p> <p data-bbox="735 797 1473 826">Magnetic homogenizer has a “M23” label engraved on the body frame.</p>
<p data-bbox="280 902 667 931">48-tube magnetic homogenizer rack</p> 	<p data-bbox="735 909 1473 992">48-tube magnetic homogenizer rack is mounted on the “M23” 48-tube magnetic homogenizer and is designed for placing and fixing reagent tubes and/or DNA samples in the wells of the rack.</p> <p data-bbox="735 992 1473 1070">The 48-tube magnetic homogenizer rack must be made of stainless steel and Trovidur EN polyvinyl chloride, in accordance with the design documentation P-240v2-01-00SB.</p> <p data-bbox="735 1070 1473 1126">The overall dimensions of the 48-tube magnetic homogenizer rack must be: 220±0,2 x 131±0,2 x 74,5±0,2 mm.</p> <p data-bbox="735 1126 1473 1155">Weight of the 48-tube magnetic homogenizer rack shall be: 900±2 g.</p> <p data-bbox="735 1155 1473 1267">Stands for a 48-tube magnetic homogenizer rack must be fixed on the “M23” homogenizer by aligning the pins on the top surface of a 48-tube magnetic homogenizer with holes in the base of a 48-tube magnetic homogenizer rack.</p> <p data-bbox="735 1267 1473 1323">The 48-tube magnetic homogenizer rack has a “T24” label engraved on its base.</p>
<p data-bbox="392 1391 555 1420">Rods dispenser</p> 	<p data-bbox="735 1397 1473 1509">Rods dispenser is installed on the worktable of the devices in versions DTstream9 L1, DTstream9 L4, DTstream12 L1, DTstream12 L4, DTstream15 L1, DTstream15 L4 and is intended for placing cartridges with magnetic “rods” in it.</p> <p data-bbox="735 1509 1473 1644">The rods dispenser consists of a base for placing the magnetic rods cartridges and magnetic tweezers for placing the rods on the test tubes. The magnetic tweezers, gripped by the nozzles of the device, capture the magnetic “rods”, carry them and place them into the test tubes with biological samples.</p> <p data-bbox="735 1644 1473 1673">The rods dispenser is made of aluminum alloy D16T.</p> <p data-bbox="735 1673 1473 1702">Dimensions of the rods dispenser: 133±0,2 x 91±0,2 x 93±0,2 mm.</p> <p data-bbox="735 1702 1473 1731">The weight of the rods dispenser unit is 770±2 g.</p> <p data-bbox="735 1731 1473 1809">The rods dispenser is fixed on the worktable of the device by means of aligning of the pins, located on the worktable of the device, and holes in the base of the rods dispenser.</p> <p data-bbox="735 1809 1473 1839">The dispenser has a “M13” label engraved on the body frame.</p>






Name	Description
<p data-bbox="331 174 513 203">1000 µl tips rack</p> 	<p data-bbox="679 181 1449 293">1000 µl tips rack is installed on the worktable of the devices in versions DTstream9 L1, DTstream9 L4, DTstream12 L1, DTstream12 L4, DTstream15 L1, DTstream15 L4 and is designed to hold platforms with 1000 µl tips when dispensing.</p> <p data-bbox="679 293 1270 322">1000 µl tips rack is made of D16T aluminum alloy and steel.</p> <p data-bbox="679 322 1359 351">Dimensions of the 1000 µl tips rack: 133±0.2 x 91±0.2 x 107±0.2 mm.</p> <p data-bbox="679 351 1062 380">Weight of 1000 µl tips rack is 715±2 g.</p> <p data-bbox="679 380 1449 517">1000 µl tips rack is placed on the worktable of the device and fixed on the table by clutching of the spring grips located in the rack body, with pins, located on the surface of the worktable. A 1000 µl tips platform is placed in the stand and held during the dispensing process. The 1000 µl tips rack has a “T11” label engraved on the body frame.</p>
<p data-bbox="338 589 507 618">200 µl tips rack</p> 	<p data-bbox="679 602 1449 714">200 µl tips rack is installed on the worktable of DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions and is designed to hold platforms with 200 µl tips when dispensing.</p> <p data-bbox="679 714 1302 743">The 200 µl tips rack is made of D16T aluminum alloy and steel.</p> <p data-bbox="679 743 1334 772">Dimensions of the 200 µl tips rack: 133±0.2 x 91±0.2 x 72±0.2 mm.</p> <p data-bbox="679 772 1088 801">The weight of 200 µl tips rack is 533±2 g.</p> <p data-bbox="679 801 1449 947">200 µl tips rack is set on the worktable of the device and fixed on the table by clutching of the spring grips located in the rack body with pins located on the surface of the worktable. A 200 µl tips platform is placed in the stand and held during the dispensing process. 200 µl tips rack has a “T12” label engraved on the body frame.</p>
<p data-bbox="248 1010 596 1070">96-well for strips, test tubes and microplates</p> 	<p data-bbox="679 1023 1449 1106">96-well for strips, test tubes and microplates is installed on the worktable of device of all versions and is designed to be placed in the wells of the rack for strips, test tubes and microplates.</p> <p data-bbox="679 1106 1442 1135">96-well for strips, test tubes and microplates is made of D16T aluminum alloy.</p> <p data-bbox="679 1135 1299 1164">Dimensions of the 96-well for strips, test tubes and microplates: 132,5±0,2 x 91±0,2 x 18±0,1 mm.</p> <p data-bbox="679 1164 1398 1193">The weight of the 96-well for strips, test tubes and microplates is 502±2 g.</p> <p data-bbox="679 1193 1449 1285">96-well for strips, test tubes and microplates has a “T01” label engraved on the body frame.</p>
<p data-bbox="242 1395 603 1456">96-well for low profile strips, test tubes and microplates</p> 	<p data-bbox="679 1408 1449 1491">96-well for low profile strips, test tubes and microplates is installed on the worktable of device of all versions and is designed to be placed in the wells of the rack for strips, test tubes and microplates.</p> <p data-bbox="679 1491 1449 1547">96-well for low profile strips, test tubes and microplates is made of D16T aluminum alloy.</p> <p data-bbox="679 1547 1449 1603">Dimensions of the 96-well for low profile strips, test tubes and microplates: 132,5±0,2 x 91±0,2 x 18±0,2 mm.</p> <p data-bbox="679 1603 1449 1659">The weight of the 96-well for low profile strips, test tubes and microplates is 414±2 g.</p> <p data-bbox="679 1659 1449 1715">96-well for low profile strips, test tubes and microplates has a “T02” label engraved on the body frame.</p>





Name	Description
<p data-bbox="316 174 630 203">4x6 rack for buffer solutions</p> 	<p data-bbox="738 188 1471 327">4x6 rack for buffer solutions is installed on the worktable of DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions and is designed for holding buffer solution containers of different volumes during the dosing process. 4x6 rack for buffer solutions is made of D16T aluminum alloy.</p> <p data-bbox="738 331 1166 360">Dimensions of the rack for buffer solutions: 132,5±0,2 x 91±0,2 x 18±0,2 mm.</p> <p data-bbox="738 392 1307 421">The weight of the 4x6 rack for buffer solutions is 504±2 g.</p> <p data-bbox="738 425 1471 504">4x6 rack for buffer solutions is installed on the worktable of the device by aligning the pins located on the worktable of the device and the holes in the base of the rack.</p> <p data-bbox="738 508 1471 562">4x6 rack for buffer solutions has a “T04” label engraved on the body frame.</p>
<p data-bbox="339 584 608 613">48-well 1.5 ml tube rack</p> 	<p data-bbox="738 598 1471 736">48-well 1.5 ml tube rack is located on the magnetic adapter M15, which is installed on the worktable of the DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions, and is designed for placing and holding test tubes with reagents and/or DNA samples in the wells.</p> <p data-bbox="738 741 1457 770">Body frame of the 48-well 1.5 ml tube rack is made of polyvinyl chloride.</p> <p data-bbox="738 775 1166 804">Dimensions of the 48-well 1.5 ml tube rack: 220±0,2 x 131±0,2 x 47±0,2 mm.</p> <p data-bbox="738 835 1267 864">The weight of the 48-well 1.5 ml tube rack is 640±2 g.</p> <p data-bbox="738 869 1471 947">48-well 1.5 ml tube rack is fixed on the magnetic adapter “M15” by aligning the pins, which are on the bottom surface of the rack, with the holes in the body of the magnetic adapter.</p> <p data-bbox="738 952 1437 981">48-well 1.5 ml tube rack has a “T17” label engraved on the body frame.</p>
<p data-bbox="280 994 668 1023">Multifunctional 48 1.5 ml tube rack</p> 	<p data-bbox="738 1008 1471 1146">Multifunctional 48 1.5 ml tube rack is located on the magnetic adapter “B 067”, which is installed on the worktable of the DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions, and is designed for placing and holding test tubes with reagents and/or DNA samples in the wells.</p> <p data-bbox="738 1151 1471 1229">Body frame of the multifunctional 48 1.5 ml tube rack is made of polyvinyl chloride.</p> <p data-bbox="738 1234 1273 1263">Dimensions of the multifunctional 48 1.5 ml tube rack: 220±0,2 x 138,5±0,2 x 50,3±0,2 mm.</p> <p data-bbox="738 1294 1374 1323">The weight of the multifunctional 48 1.5 ml tube rack is 710±2 g.</p> <p data-bbox="738 1328 1471 1406">Multifunctional 48 1.5 ml tube rack is fixed on the magnetic adapter “M15” by aligning the pins, which are on the bottom surface of the rack, with the holes in the body of the magnetic adapter.</p> <p data-bbox="738 1411 1471 1464">Multifunctional 48 1.5 ml tube rack has a “T16” label engraved on the body frame.</p>
<p data-bbox="347 1458 600 1523">Magnetic adapter for 48 1.5 ml tube rack</p> 	<p data-bbox="738 1471 1471 1610">Magnetic adapter for 48 1.5 ml tube rack is installed on the worktable of the DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions and is designed to hold a “T17” rack for 48 1.5 ml tubes.</p> <p data-bbox="738 1615 1471 1693">Body frame of the magnetic adapter for 48 1.5 ml tube rack is made of polyvinyl chloride.</p> <p data-bbox="738 1697 1471 1776">The magnetic adapter for 48 1.5 ml tube rack has 48 cylindrical permanent magnets. The magnetic field of the magnets interacts with magnetic particles and isolated DNA molecules in the tubes, depositing them on the walls of the tubes.</p> <p data-bbox="738 1780 1326 1809">Dimensions of the magnetic adapter for 48 1.5 ml tube rack: 220±0,2 x 131±0,2 x 34±0,2 mm.</p> <p data-bbox="738 1841 1437 1870">The weight of the magnetic adapter for 48 1.5 ml tube rack is 1046±2 g.</p> <p data-bbox="738 1874 1471 1953">Magnetic adapter for 48 1.5 ml tube rack is fixed on the worktable by aligning the pins, which are on the bottom surface of the rack, with the holes in the body of the magnetic adapter.</p> <p data-bbox="738 1957 1471 2011">Magnetic adapter for 48 1.5 ml tube rack has a “M15” label engraved on the body frame.</p>





Name	Description
<p data-bbox="236 174 609 235">Adapter with light pointer for tube arrangement</p> 	<p data-bbox="676 181 1445 293">Adapter with light pointer for tube arrangement is installed on the worktable of the DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions and is designed for placement and retention of 48 1.5 ml tubes.</p> <p data-bbox="676 293 1445 344">Body frame of the adapter with light pointer for tube arrangement is made of polyvinyl chloride.</p> <p data-bbox="676 344 1321 374">Dimensions of the adapter with light pointer for tube arrangement:</p> <p data-bbox="676 374 1002 403">228±0,2 x 139±0,2 x 47±0,2 mm.</p> <p data-bbox="676 403 1422 432">The weight of the adapter with light pointer for tube arrangement is 930±2 g.</p> <p data-bbox="676 432 1445 591">Under each well of the adapter with light pointer for tube arrangement there is a light-emitting diode, which illumination should indicate to the operator the location of the tube. The LEDs are powered and the control signal from the device control software is supplied through a cable that connects the USB connector on the device panel to the micro-USB connector in the adapter body frame.</p> <p data-bbox="676 591 1445 674">Adapter with light pointer for tube arrangement is fixed on the worktable by aligning the pins located on the worktable of the device and the holes in the base of the device.</p> <p data-bbox="676 674 1445 725">Adapter with light pointer for tube arrangement has a “P-154” label engraved on the body frame.</p>
<p data-bbox="261 775 584 804">Adapter for reagent cartridges</p> 	<p data-bbox="676 781 1445 893">Adapter for reagent cartridges is installed on the worktable of the devices in versions DTstream9 L1, DTstream9 L4, DTstream12 L1, DTstream12 L4, DTstream15 L1, DTstream15 L4 and is designed for placing two reagent cartridges in it.</p> <p data-bbox="676 893 1302 922">Adapter for reagent cartridges is made of D16T aluminum alloy.</p> <p data-bbox="676 922 1152 952">Dimensions of the adapter for reagent cartridges:</p> <p data-bbox="676 952 991 981">133±0,2 x 92±0,2 x 14±0,2 mm.</p> <p data-bbox="676 981 1251 1010">The weight of the adapter for reagent cartridges is 274±2 g.</p> <p data-bbox="676 1010 1445 1093">Adapter for reagent cartridges is fixed on the worktable by aligning the pins located on the worktable of the device and the holes in the base of the device.</p> <p data-bbox="676 1093 1445 1144">Adapter for reagent cartridges has a “M14” label engraved on the body frame.</p>
<p data-bbox="271 1171 574 1200">384-well microplate adapter</p> 	<p data-bbox="676 1178 1445 1261">384-well microplate adapter is installed on the worktable of the DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions and is designed to fix 384-well microplate in it.</p> <p data-bbox="676 1261 1425 1290">384-well microplate adapter is made of D16T aluminum alloy and 65G steel.</p> <p data-bbox="676 1290 1142 1319">Dimensions of the 384-well microplate adapter:</p> <p data-bbox="676 1319 1043 1348">133,5±0,2 x 91,2±0,2 x 12,5±0,2 mm.</p> <p data-bbox="676 1348 1240 1377">The weight of the 384-well microplate adapter is 195±2 g.</p> <p data-bbox="676 1377 1445 1460">384-well microplate adapter is fixed on the worktable by aligning the pins located on the worktable of the device and the holes in the base of the adapter.</p> <p data-bbox="676 1460 1410 1489">384-well microplate adapter has a “T03” label engraved on the body frame.</p>
<p data-bbox="268 1588 577 1617">192-well microplates adapter</p> 	<p data-bbox="676 1594 1445 1677">192-well microplates adapter is installed on the worktable of the DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions and is designed to fix 192-well microplate in it.</p> <p data-bbox="676 1677 1382 1706">192-well microplates adapter is made of organic sheet glass (Plexiglass).</p> <p data-bbox="676 1706 1150 1736">Dimensions of the 192-well microplates adapter:</p> <p data-bbox="676 1736 991 1765">133±0,2 x 86±0,2 x 15±0,2 mm.</p> <p data-bbox="676 1765 1238 1794">The weight of the 192-well microplates adapter is 63±2 g.</p> <p data-bbox="676 1794 1445 1877">192-well microplates adapter is fixed on the worktable by aligning the pins located on the worktable of the device and the holes in the base of the adapter.</p> <p data-bbox="676 1877 1422 1906">192-well microplates adapter has a “T08” label engraved on the body frame.</p>



Name	Description
<p data-bbox="323 174 627 203">192-well microplate adapter</p> 	<p data-bbox="735 181 1471 293">192-well microplate adapter is installed on the worktable of the DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4 versions and is designed to fix 192-well microplate in it.</p> <p data-bbox="735 293 1422 322">192-well microplate adapter is made of Trovidur ER polymer material.</p> <p data-bbox="735 322 1198 351">Dimensions of the 192-well microplate adapter: 132,5±0,2 x 91±0,2 x 18,1±0,2 mm.</p> <p data-bbox="735 374 1299 403">The weight of the 192-well microplate adapter is 107±2 g.</p> <p data-bbox="735 403 1471 483">192-well microplate adapter s fixed on the worktable by aligning the pins located on the worktable of the device and the holes in the base of the adapter.</p> <p data-bbox="735 483 1466 512">192-well microplate adapter has a “T07” label engraved on the body frame.</p>
<p data-bbox="280 629 670 683">Adapter for used material collection container</p> 	<p data-bbox="735 631 1471 712">Adapter for used material collection container is installed on the worktable of all versions of the device and is designed to place and hold the MK-01 container for waste material collection.</p> <p data-bbox="735 712 1471 768">Adapter for used material collection container is made of D16T aluminum alloy.</p> <p data-bbox="735 768 1362 797">Dimensions of the adapter for used material collection container: 143,9±0,2 x 108,4±0,2 x 13±0,2 mm.</p> <p data-bbox="735 819 1466 848">The weight of the adapter for used material collection container is 184±2 g.</p> <p data-bbox="735 848 1471 929">Adapter for used material collection container s fixed on the worktable by aligning the pins located on the worktable of the device and the holes in the base of the device.</p> <p data-bbox="735 929 1471 987">Adapter for used material collection container has a “K05” label engraved on the body frame.</p>
<p data-bbox="280 1126 670 1180">Stand for DTstream liquid handling station</p> 	<p data-bbox="735 1128 1471 1187">Stand for DTstream liquid handling station is designed for the installation of a station on its working surface.</p> <p data-bbox="735 1187 1471 1292">The countertops of both models are made of Trovidur EN polyvinyl chloride (plastic sheet) and have an opening on the surface to allow the discharge of waste materials (tips) from the device into the waste material container inside the stand.</p> <p data-bbox="735 1292 1471 1373">When mounting the device on a stand, the outlet for discharging the station waste materials (tips) must be aligned with the intake opening on the stand table top.</p> <p data-bbox="735 1373 1471 1431">If necessary Inside the stand under the intake for discharge of waste materials (tips), install the tips receptacle “K12-12-12-00-00”.</p> <p data-bbox="735 1431 1471 1541">Structurally the stand is made as an all-welded frame from a steel pipe of rectangular section, non-working surfaces are made of metal sheet, the swing doors are on hinges, the drawer in the stand “CtDy9v2-00-00-00-00” is moved along the rails.</p> <p data-bbox="735 1541 1471 1599">The feet of the stands provide height adjustment to ensure the horizontal position of the stand countertop during operation of the device.</p> <p data-bbox="735 1599 1016 1628">The dimensions of the stand: 1200±0,2 x 745±0,2 x 900±0,2 mm.</p> <p data-bbox="735 1650 1107 1680">The weight of the stand is 103±0,5 kg.</p> <p data-bbox="735 1680 1198 1709">The load capacity of the stand is at least 250 kg.</p> <p data-bbox="735 1709 1294 1738">On the body of the stand is a nameplate with the marking.</p>

Name	Description
<p data-bbox="316 174 536 203">Waste tip receptacle</p> 	<p data-bbox="683 174 1445 259">Waste tip receptacle is installed inside the stand under the device “CtDy9v2-00-00-00-00” and is designed to place a container on it for receiving waste materials.</p> <p data-bbox="683 259 1445 315">Waste tip receptacle is made of D16T aluminum alloy, Trovidur EN polyvinyl chloride, caprolon.</p> <p data-bbox="683 315 1406 344">Dimensions of the waste tip receptacle: $441\pm 0,2 \times 220\pm 0,2 \times 327\pm 0,2$ mm.</p> <p data-bbox="683 344 1166 374">The weight of the waste tip receptacle: $475\pm 0,5$ g.</p> <p data-bbox="683 374 1310 403">The body frame of the device has a “K12-12-12-00-00” marking.</p>
<p data-bbox="276 589 576 651">Protective cap for magnetic tweezers</p> <p data-bbox="304 667 547 696">Cat. No. C-DTS-C002</p> 	<p data-bbox="683 600 1445 685">Protective cap for magnetic tweezers is installed on the rods dispenser “M13” and is intended for single use in order to protect the working surface of the magnetic tweezers before each dispensing.</p> <p data-bbox="683 685 1414 714">Protective cap for magnetic tweezers is made of transparent polypropylene.</p> <p data-bbox="683 714 1310 743">Dimensions of the protective cap: $61\pm 0,2 \times 7\pm 0,2 \times 6,5\pm 0,2$ mm.</p> <p data-bbox="683 743 1102 772">The weight of the protective cap is $1\pm 0,2$ g.</p>
<p data-bbox="268 1061 584 1124">Replaceable funnel for waste materials</p> <p data-bbox="304 1140 547 1169">Cat. No. C-DTS-C001</p> 	<p data-bbox="683 1084 1445 1169">Replaceable funnel for waste materials is installed on the devices of all versions and is designed for single use for the safe discharge of waste materials (tips).</p> <p data-bbox="683 1169 1445 1225">Replaceable funnel for waste materials is made of transparent plastic PET (polyethylene terephthalate) with a thickness of 0.5 mm.</p> <p data-bbox="683 1225 1326 1254">Dimensions of the funnel: $220\pm 0,2 \times \text{Ø}54,5\pm 0,2 \times \text{Ø}38,6\pm 0,2$ mm.</p> <p data-bbox="683 1254 1334 1283">The weight of the replaceable funnel for waste materials: $12\pm 0,5$ g.</p> <p data-bbox="683 1283 1445 1368">Replaceable funnel for waste materials is installed in the receiving socket of the liquid handling station stand “CtDy9v2-00-00-00-00” and sends the waste material to the waste materials reception container.</p>
<p data-bbox="316 1520 536 1550">CAN interface cable</p> 	<p data-bbox="683 1532 1445 1727">CAN interface cable is designed for DTstream9 L1, DTstream9 L4, DTstream12 L1, DTstream12 L4, DTstream15 L1, DTstream15 L4 version with 48-tube magnetic homogenizer. The cable consists of connector type MINI DIN-8, straight; connector type “C091-U 250V/5A,8Pol Male T 3504 005 U”; connection cable in polyurethane isolation. Cable dimensions: $450\pm 10 \times \text{d}6$ mm. The cable connects 48-tube magnetic homogenizer to the CON connector on the switching panel.</p>

Name	Description
<p data-bbox="331 174 616 203">Power cable (three-prong)</p> 	<p data-bbox="740 181 1465 237">Power cable (three-wire) is designed to supply power from the mains to the product.</p> <p data-bbox="740 239 884 264">Specifications:</p> <p data-bbox="740 266 1050 291">Socket: molded PVS 45P, black</p> <p data-bbox="740 293 858 318">16 A, 250 V</p> <p data-bbox="740 320 1031 344">Plug: Molded PVS 45P, black</p> <p data-bbox="740 347 935 371">Terminal: 2x24 mm</p> <p data-bbox="740 374 1238 398">Wire: H05 W-F 3G 0.75 mm² GTSA-3, OD6.8 mm</p> <p data-bbox="740 400 919 425">Length: 1800 mm.</p>
<p data-bbox="309 472 639 501">Ethernet communication cable</p> 	<p data-bbox="740 479 1465 535">Ethernet communication cable is designed for communication with a personal computer</p> <p data-bbox="740 537 884 562">Specifications:</p> <p data-bbox="740 564 1453 620">Patch Cord UTP 5e cat. "HK-SC5EUTP-RD-2.0" with RJ-45 connectors, 26AWG/0,4 mm</p> <p data-bbox="740 622 1007 647">Cable: cat5e UTP 26AWG.</p> <p data-bbox="740 649 1091 674">PVC jacket, cable diameter: 5.5 mm</p> <p data-bbox="740 676 839 701">Cap: PVC</p> <p data-bbox="740 703 1118 728">RJ-45 plug: polycarbonate, cat5e, 8p8c</p> <p data-bbox="740 730 863 754">Category: 5e</p> <p data-bbox="740 757 983 781">Design: unshielded, UTP</p> <p data-bbox="740 784 1027 808">Connector format: RJ45/8p8c</p> <p data-bbox="740 810 975 835">Cable diameter: 5.5 mm</p> <p data-bbox="740 837 884 862">Length: 2.0 m.</p>
<p data-bbox="316 898 632 927">Fuses (10A, 5x20mm, 250V)</p> 	<p data-bbox="740 904 1465 960">Fuses (10A, 5x20mm, 250V) are designed to protect the electronic unit against overloading.</p> <p data-bbox="740 963 884 987">Specifications:</p> <p data-bbox="740 990 1038 1014">Fuse type: cylindrical, ceramic</p> <p data-bbox="740 1016 932 1041">Rated current: 10 A</p> <p data-bbox="740 1043 979 1068">Rated voltage: AC 250V</p> <p data-bbox="740 1070 959 1095">Dimensions: 5x20 mm</p> <p data-bbox="740 1097 1190 1122">Maximum breaking capacity: 1500 AC 250 V.</p>
<p data-bbox="336 1173 608 1252">1000 µl filter tips Cat. No. C-DTS-T1000F</p> 	<p data-bbox="740 1180 1465 1263">1000 µl filter tips are designed to perform dosing of solutions and reagents for DTstream9 L1, DTstream9 L4, DTstream12 L1, DTstream12 L4, DTstream15 L1, DTstream15 L4.</p> <p data-bbox="740 1265 1465 1321">The tips are placed in the platform in the amount of 96 pieces. A filter is placed inside the tip.</p> <p data-bbox="740 1323 1422 1348">Dimensions of the platform for tips: 135±0,2 x 97±0,2 x 105±0,2 mm.</p> <p data-bbox="740 1350 1465 1406">The platform with tips is installed in the rack for tips 1000 µl "T11", placed on the table of the liquid handling station.</p>
<p data-bbox="344 1554 600 1632">1000 µl tips Cat. No. C-DTS-T1000</p> 	<p data-bbox="740 1561 1465 1644">1000 µl tips are designed for dosing solutions and reagents for devices in versions DTstream9 L1, DTstream9 L4, DTstream12 L1, DTstream12 L4, DTstream15 L1, DTstream15 L4.</p> <p data-bbox="740 1646 1342 1671">The tips are placed in the platform in the amount of 96 pieces.</p> <p data-bbox="740 1673 1422 1697">Dimensions of the platform for tips: 135±0,2 x 97±0,2 x 105±0,2 mm.</p> <p data-bbox="740 1700 1465 1756">The platform with the tips is installed in the rack for tips "T11", placed on the table of the liquid handling station.</p>

Name	Description
<p>200 µl filter tips Cat. No. C-DTS-T200F</p> 	<p>200 µl filter tips are designed for dosing of solutions and reagents for DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream12 M4, DTstream15 M1, DTstream15 M4.</p> <p>The tips are placed on the platform in the amount of 96 pieces. Inside the tip there is a filter.</p> <p>Dimensions of the platform for tips: 136±0,2 x 98±0,2 x 70±0,2 mm.</p> <p>The platform with tips is installed in the rack for tips “T12”, placed on the worktable of the device.</p>
<p>200 µl tips Cat. No. C-DTS-T200</p> 	<p>200 µl tips are designed for dosing solutions and reagents for DTstream9 M1, DTstream9 M4, DTstream12 M1, DTstream9 M4, DTstream15 M1, DTstream15 M4.</p> <p>The tips are placed in the platform in the amount of 96 pieces.</p> <p>Dimensions of the platform for tips: 136±0,2 x 98±0,2 x 70±0,2 mm.</p> <p>The platform with tips is installed in the rack for tips “T12”, placed on the worktable of the device.</p>
<p>96-well microplate Cat. No. C-DTS-P096/01</p> 	<p>96-well microplate is designed to be filled with liquid reagents and samples of biological material for further use in PCR amplification.</p> <p>Dimensions according to specification:</p> <p>Width: 125.11±0.25 mm</p> <p>Depth: 83.22±0.25 mm</p> <p>Height: 20.80±0.25 mm</p> <p>Well depth: 20.20±0.15 mm</p> <p>Diameter of the wells: 5.46±0.10 mm</p> <p>The microplate is placed on the worktable of the device in the “T01” rack.</p>
<p>96-well semi-skirted microplate Cat. No. C-DTS- P096/02</p> 	<p>96-well semi-skirted microplate is designed to be filled with liquid reagents and samples of biological material for further use in PCR amplification.</p> <p>Dimensions according to manufacturer's specifications:</p> <p>Width: 124.26±0.25 mm</p> <p>Depth: 83.97±0.25 mm</p> <p>Height: 20.70±0.25 mm</p> <p>Well depth: 20.20±0.10 mm</p> <p>Well diameter: 5.46±0.10 mm</p> <p>The microplate is placed on the worktable of the device in the “T01” rack.</p>

Name	Description
<p>192-well microplate Cat. No. C-DTS-P192/01</p> 	<p>192-well microplate is designed to be filled with liquid reagents and samples of biological material for further use in PCR amplification. Dimensions according to specification: Width: 80.0(-0.25)mm Depth: 60.0(-0.25)mm Height: 10.10±0.25 mm Well depth: 9.60±0.10 mm Well diameter: 3.10±0.10 mm The microplate is placed on the worktable of the device in the “T08” adapter.</p>
<p>192-well semi-skirted microplate Cat. No. C-DTS-P192/02</p> 	<p>192-well semi-skirted microplate is designed to be filled with liquid reagents and samples of biological material for further use in PCR amplification. Dimensions according to specifications: Width: 80.0(-0.25)mm Depth: 60.0(-0.25)mm Height: 10.10±0.25 mm Well depth: 9.60±0.10 mm Well diameter: 3.10±0.10 mm The microplate is placed on the worktable of the device in “T07”, “T08” adapters.</p>
<p>384-well microplate Cat. No. C-DTS-P384</p> 	<p>384-well microplate is designed to be filled with liquid reagents and samples of biological material for further use in PCR amplification. Dimensions according to specification: Width: 127.76±0.25mm. Depth: 85.48±0.25 mm Height: 10,60±0,25 mm Well depth: 9.60±0.10 mm Well diameter: 3.10±0.10 mm The microplate is placed on the worktable of the device in the “T03” adapter.</p>
<p>Container for waste material collection on the device worktable</p> 	<p>As a container for waste material collection it is recommended to use a medical product – “Containers for medical waste and consumables (container for collection, storage, transportation and disposal of sharps waste (yellow) with a volume of 1.3 l)”. The container for waste material collection is installed on the “K05” adapter, placed on the worktable of the device.</p>

Name	Description
<p data-bbox="304 176 587 210">Microplate seal (package)</p> <p data-bbox="325 226 566 255">Cat. No. C-DTS-F001</p>  <p data-bbox="432 562 459 584">or</p> <p data-bbox="330 607 561 636">Microplate seal (roll)</p> 	<p data-bbox="716 185 1449 241">Microplate seal is designed to thermally seal microplate wells after filling them with reagent solutions and/or DNA samples.</p> <p data-bbox="716 241 1310 271">The microplate seal can be supplied in cut sheets or on a roll.</p> <p data-bbox="716 271 1374 300">A package of 100 sheets of seal is packed in a polyethylene zip bag.</p> <p data-bbox="716 300 1102 329">Dimensions of seal sheets: 125x78 mm.</p> <p data-bbox="716 329 1150 358">Dimensions of the seal roll: 610 m x 78 mm.</p>

Attention! Components from other manufacturers must not be used.

The manufacturer is not responsible for the performance of the DTstream liquid handling station if the customer uses third-party components.