



For research use only

## HIV PCR detection Kit USER MANUAL



"DNA-Technology, Research & Production" LLC

Russia, 142281, Moscow Region,

Protvino, 20 Zheleznodorozhnaya Street,

Phone/fax: +7(495) 640.17.71

+7(4967)31.06.70,

E-mail: [protvino@dna-technology.ru](mailto:protvino@dna-technology.ru), [mail@dna-technology.ru](mailto:mail@dna-technology.ru)

Customer service department: +7(495)640.16.93

E-mail: [hotline@dna-technology.ru](mailto:hotline@dna-technology.ru),

[www.dna-technology.ru](http://www.dna-technology.ru)



R3-P609-23/9EU  
R3-P609-S3/9EU  
R3-P609-24/9EU  
F3-P609-21/1EU  
F3-P609-51/1EU  
F3-P609-22/1EU  
F3-P609-52/1EU



311-3.2024.04.22

## Table of contents

<b>1. INTENDED USE</b>	<b>4</b>
<b>2. METHOD</b>	<b>4</b>
<b>3. CONTENT</b>	<b>5</b>
<b>4. REAGENTS AND EQUIPMENT REQUIRED BUT NOT PROVIDED</b>	<b>6</b>
<b>5. WARNINGS AND PRECAUTIONS</b>	<b>7</b>
<b>6. RNA EXTRACTION PROTOCOL</b>	<b>9</b>
<b>7. REVERSE TRANSCRIPTION PROTOCOL</b>	<b>11</b>
<b>8. PCR PROTOCOL</b>	<b>12</b>
<b>9. CONTROLS</b>	<b>15</b>
<b>10. DATA ANALYSIS</b>	<b>12</b>
<b>11. TROUBLESHOOTING</b>	<b>16</b>
<b>12. STORAGE AND HANDLING REQUIREMENTS</b>	<b>17</b>
<b>13. SPECIFICATIONS</b>	<b>17</b>
<b>14. QUALITY CONTROL</b>	<b>18</b>
<b>15. KEY TO SYMBOLS</b>	<b>18</b>
<b>Annex A</b>	<b>20</b>

## 1. INTENDED USE

The **HIV PCR detection Kit** is intended for research and diagnostic applications. The **HIV PCR detection Kit** is an *in vitro* Nucleic Acid Test (NAT) based pathogen detection product. The **HIV PCR detection Kit** is designed to detect Human Immunodeficiency Virus type 1 (HIV-1) nucleic acids in human blood plasma.

The **HIV PCR detection Kit** can be used in clinical practice for HIV diagnostics.

## 2. METHOD

The implemented method is based on viral RNA reverse transcription followed by PCR amplification of the obtained cDNA.

The detection can be performed in each of two variants: real-time, and end-point fluorescent detection (FLASH technology).

The **HIV REAL-TIME PCR Detection Kit** is based on fluorescent modification of the PCR method. The PCR-mix contains target-specific probes bearing reporter and quencher molecules. Once hybridized to a target sequence, the probes become activated. As a result of activation fluorescence increases proportionally to target sequence amplification. The intensity of fluorescence is measured at every cycle of reaction with a Real-time PCR thermal cycler data collection unit and analyzed with the software provided. The **HIV FLASH PCR Detection Kit** is based on the same principle but the fluorescence is measured only once after reaction.

The automatic analysis for **HIV REAL-TIME PCR Detection Kit** is available on “DNA-Technology” made instruments: DTlite and DTprime Thermal Cyclers (see the catalogue at [www.dna-technology.ru/en](http://www.dna-technology.ru/en) to see available supply options) and for **HIV FLASH PCR Detection Kit** on Gene or Gene-4 Fluorescence Readers

**REF** O-GENE-EU, O-GENE4-EU.

The **HIV REAL-TIME PCR Detection Kit** **REF** R3-P609-23/9EU, R3-P609-S3/9EU and R3-P609-24/9EU is also approved for use with iQ (Bio-Rad Laboratories) and Rotor-Gene Q (Qiagen) real-time thermal cyclers.

The Kit can be supplied in either separate (1x96) or stripped (8x12) tubes (**REF** R3-P609-23/9EU, R3-P609-24/9EU and R3-P609-S3/9EU respectively).

### 3. CONTENT

The Kit consists of three parts, namely: PREP-NA DNA/RNA Extraction Kit, Reverse Transcription Kit and HIV PCR detection Kit (either Real-time or FLASH variant) each is providing the crucial step of the assay.

Table 1. PREP-NA DNA/RNA Extraction Kit<sup>1</sup>

Reagent	Description	Total volume	Amount
Lysis buffer	Light blue slightly foaming liquid	30 mL	1 vial
Precipitation buffer	Colorless liquid	40 mL	1 vial
Washout solution 1	Colorless liquid	50 mL	1 vial
Washout solution 2	Colorless liquid	30 mL	1 vial
Dilution buffer	Colorless liquid	5 mL (1.25 mL in each tube)	4 tubes
Negative control (C-)	Colorless liquid	3 mL (1.5 mL in each tube)	2 tubes
Internal control (RNA-IC)	Colorless liquid	1 mL	1 tube

Table 2. Reverse Transcription Kit

Reagent	Description	Total volume	Amount
RT-buffer	Colorless liquid	200 µL	1 tube
RT-HAV+HCV+HDV+HGV+HIV primers and dNTP's	Colorless liquid	100 µL	1 tube
Reverse transcriptase	Colorless viscous liquid	50 µL	1 tubes

---

<sup>1</sup> - can be included into the kit if requested.

Table 3. HIV PCR detection Kit

Reagent	Description	Total volume	Amount
Paraffin sealed PCR-mix	Colorless liquid and white waxy fractions	1.92 mL or 2,0 mL (0.02 mL per tube)	96 or 100 separate or stripped tubes of 0,2 or 0,5 mL
TECHNO Taq-polymerase	Colorless viscous liquid	50 µL	1 tube
PCR-buffer	Colorless liquid	1 mL (0.5 mL in each tube)	2 tubes
Positive control (C+)	Colorless liquid	150 µL	1 tube
Mineral oil (not supplied in Kit for Rotor-Gene Q)	Colorless viscous liquid	2 mL (1 mL in each tube)	2 tubes

The approximate total time needed to perform the assay is 4 hours.

The *PREP-NA* DNA/RNA Extraction Kit is sufficient for extraction of 100 samples.

The **HIV PCR detection Kit** sufficient to test 96/100 samples including IC, negative and positive control samples.

#### 4. REAGENTS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

##### 4.1 Specimen collection

The blood samples should be collected in 2.0 or 4.0 mL Vacuette type tubes with EDTA in 2.0 mg/mL final concentration. The sodium citrate anticoagulant is also applicable.



The use of heparin anticoagulant is not allowed.

##### 4.2 RNA extraction and PCR

Class II biological safety cabinet;

UV PCR cabinet;

Vortex mixer;

Aspirator with trap flask to remove supernatants;

0.2, 0.5 and 1.5 mL tubes;

PCR tube rack for 0.2, 0.5 and 1.5 mL tubes;

Single channel pipettes (volume range 2-20 µL, 20-200 µL, 200-1000 µL);

RNase and DNase free filtered pipette tips (volume range 20 µL, 200 µL, 1000 µL);

Pipette tip without filter (for electrical laboratory aspirator);

Powder-free surgical gloves;

Disinfectant solution;

Container for used pipette tips;

High speed centrifuge (RCF 16 000 x g);

Thermostat (temperature range from 40 °C to 95 °C);

Real-time PCR thermal cycler (for **HIV REAL-TIME PCR Detection Kit**); if using “DNA-Technology” made

DTprime **REF** O-DTPRIME4M1-EU, O-DTPRIME5M1-EU, DTlite **REF** O-DTLITE4S1-EU, O-DTLITE5S1-EU, software version must be not lower than 7.3; «HIV.ini» file<sup>2</sup>.

Tercyc Conventional PCR Thermal Cycler (**REF** O-TP4-EU) or equivalent (for **HIV FLASH PCR Detection Kit**);

Gene or Gene-4 Fluorescence Reader (**REF** O-GENE-EU, O-GENE4-EU) or Ala1/4 fluorescence reader or equivalent (for **HIV FLASH PCR Detection Kit**); for Gene Fluorescence Reader software version should be not lower than 3.3, for Gene-4 Fluorescence Reader software version must be not lower than 4.4.0.8.

## 5. WARNINGS AND PRECAUTIONS

As part of industrial and personal hygiene and general safety practices, avoid all unnecessary exposure to the chemical components of this kit and ensure prompt removal from skin, eyes, and clothing upon contact.

Significant health effects are **NOT** anticipated from routine use of this kit when adhering to the instructions listed in the current manual.

### 5.1 General warnings and precautions

Handle and dispose all biological samples as if they were able to transmit infective agents. Avoid direct contact with the biological samples. Avoid producing spills or aerosol. Any material coming in contact with the biological samples must be treated for at least 30 minutes with disinfecting solution or autoclaved for 1 hour at 121°C before disposal.

Handle and dispose all reagents and all materials used to carry out the assay as if they were able to transmit infective agents. Avoid direct contact with the reagents. Avoid producing spills or aerosol. Waste

---

<sup>2</sup> - - the latest recommended version can be download from the company's web-site <http://www.dna-technology.ru/eng/support/>

must be handled and disposed according to adequate safety measures. Disposable combustible material must be incinerated. Liquid waste containing acids or bases must be neutralized before disposal.

Wear suitable protective clothes and gloves and protect eyes and face.

Never pipette solutions by mouth.

Do not eat, drink, smoke or apply cosmetic products in the work areas.

Carefully wash hands after handling samples and reagents.

Dispose of leftover reagents and waste in compliance with the regulations in force.

Read the **User manual** provided with the kit before running the assay.

While running the assay follow the instructions listed in the **User manual**.

Do not use the kit after the expiry date provided.

Only use the reagents provided in the kit and those recommended by manufacturer.

Do not mix reagents from different batches.

Do not use reagents from third party manufacturers' kits.

## **5.2 Warnings and precautions for molecular biology**

Molecular biology procedures, such as nucleic acids extraction, reverse transcription, amplification and detection require qualified staff to avoid the risk of erroneous results, especially due to the degradation of nucleic acids contained in the samples or sample contamination by amplification products.

Equip separate areas for the extraction/preparation of amplification reactions and for the amplification/detection of amplification products. Never introduce an amplification product in the area designed for extraction/preparation of amplification reactions.

Wear lab coats, gloves and tools, which are exclusively employed for the extraction/preparation of the amplification reaction and for the amplification/detection of the amplification products. Never transfer lab coats, gloves and tools from the area designed for amplification/detection of the amplification products to the area designed for extraction/preparation of amplification reactions.

The samples must be exclusively employed for certain type of analysis. Samples must be handled under a laminar flow hood. Tubes containing different samples must never be opened at the same time. Pipettes used to handle samples must be exclusively employed for this specific purpose. The pipettes must be of the positive dispensation type or be used with aerosol filter tips. The tips employed must be sterile, free from the DNases and RNases, free from DNA and RNA.

The reagents must be handled under a laminar flow hood. The reagents required for amplification must be prepared in such a way that they can be used in a single session. Pipettes used to handle reagents must be exclusively employed for this specific purpose. The pipettes must be of the positive dispensation type or be used with aerosol filter tips. The tips employed must be sterile, free from the DNases and RNases, free from DNA and RNA.

Amplification products must be handled in such a way as to reduce dispersion into the environment as much as possible, in order to avoid the possibility of contamination. Pipettes used to handle amplification products must be exclusively employed for this specific purpose.

### 5.3 Emergency actions

**Inhalation:** Inhalation of the Master Mix contained within this kit is unlikely, however care should be taken.

**Eye Contact:** If any component of this kit enters the eyes, wash eyes gently under potable running water for 15 minutes or longer, making sure that the eyelids are held open. If pain or irritation occurs, obtain medical attention.

**Skin Contact:** If any component of this kit contacts the skin and causes discomfort, remove any contaminated clothing. Wash affected area with plenty of soap and water. If pain or irritation occurs, obtain medical attention.

**Ingestion:** If any component of this kit is ingested, wash mouth out with water. If irritation or discomfort occurs, obtain medical attention.

## 6. RNA EXTRACTION PROTOCOL

The **HIV PCR detection Kit** is designed to detect cDNA extracted from blood plasma. Shake the tube containing blood sample thoroughly to mix the blood and anticoagulant.



Using of heparin as anticoagulant is not allowed.



The overall storage of the sample should not exceed 6 hours.

The transportation and storage temperature from collecting the sample till analysis should be between 2 °C and 8 °C.



Whole blood cannot be frozen.

- 6.1** To obtain the plasma spin the tubes with blood at 800-1600 x g (corresponds to 3000 rpm on Eppendorf Centrifuge 5424) for 20 min at room temperature (between 18 °C and 25 °C).



Relative centrifugal force (RCF or g) depends on rotation frequency and centrifugation radius (Annex A). To establish if your centrifuge meets the requirements apply to the exploitation manual for centrifuge.

- 6.2** Take the upper fraction (plasma) with an automatic sampler and put it into the new 1.5 mL tube. The blood plasma can be stored at minus20 °C for 3 months.





The lysis buffer can contain the precipitate. Keep it at 65 °C until precipitate fully dissolved.



At this step of assay use only disposable pipette tips with filter and are RNase and DNase free.



If the plasma was frozen, prior to extraction thaw it at room temperature or at temperature from 2 °C to 8 °C. Vortex the thawed plasma and centrifugate at 16000 x g for 5 min. Take the volume needed for RNA extraction without touching the precipitate. Not more than 3 freezing/thawing cycles are allowed.



To increase the sensitivity of analysis it is necessary to carry out preliminary ultracentrifugation at 16500 rpm (24000 x g) for 1 hour at temperature between 2 °C and 8 °C.

**6.3** Mark the required number of 1.5 mL tubes by the following scheme: for each test sample and for negative control (“C-”).

For example: if you need to test 10 samples, mark 11 tubes (10 for the samples, 1 for C-).

**6.4** Add 10 µL of the premixed internal control (RNA-IC) in each tube.

**6.5** Add 300 µL of the lysis buffer avoiding contact of the pipette tip with an edge of the tube. Close the tubes.



Open the tube, add sample, then close the tube before proceeding to the next DNA sample to prevent contamination.

**6.6** Add 100 µL of the premixed blood plasma sample into the marked tubes. Do not add samples to the “C-” tube.

**6.7** Add 100 µL of the “C-” into corresponding tube.

**6.8** Close the tubes and vortex them maximally intensively for 3–5 sec twice.

**6.9** Incubate the tubes for 15 min at 65 °C, spin down the drops at 16000 x g for 30 s at room temperature.

**6.10** Add 400 µL of the precipitation buffer into each tube. Close the tubes and vortex them maximally intensively for 3-5 sec twice.

**6.11** Spin the tubes at 16000 x g for 15 min at room temperature.

**6.12** Remove the supernatant avoiding contact of the pipette tip with the precipitate. Use new tip for each sample.

**6.13** Add 500 µL of the washout solution №1 to the precipitate. Vortex the tubes for 3–5 sec. Invert the tube up and down washing of the cap of the tube. It is necessary to carry out this procedure with each tube individually

**6.14** Spin the tubes at 16000 x g for 5 min at room temperature.

**6.15** Remove the supernatant avoiding contact of the pipette tip with the precipitate. Use new tip for each sample.

- 6.16** Add 300  $\mu\text{L}$  of the washout solution №2 to the precipitate. Vortex the tubes for 3–5 sec. Gently invert the tube up and down washing of the walls and the cap of the tube. It is necessary to carry out this procedure with each tube individually.
- 6.17** Spin the tubes at 16000 x g for 5 min at room temperature.
- 6.18** Remove the supernatant avoiding contact of the pipette tip with the precipitate. Use new tip for each sample.
- 6.19** Open the tubes and dry the precipitate at 65 °C for 5 min strictly.
- 6.20** Add 16.5  $\mu\text{L}$  of the dilution buffer to the precipitate. Spin down the drops for 3–5 sec.
- 6.21** Incubate the tubes for 10 min at 65 °C, vortex them maximally intensively.
- 6.22** Spin down the drops at 16000 x g for 30 sec.



The RNA preparation is ready. RNA should be use immediately for reverse transcription reaction, RNA sample shouldn't be stored!



For several infections testing add 50  $\mu\text{L}$  of the dilution buffer to the precipitate. Keep in mind that increasing of the buffer volume results in the proportional dilution of the sample and reducing analysis sensitivity. Analytical sensitivity 200 copies/mL is verified for dilution buffer volume equal 16.5  $\mu\text{L}$ .

## 7. REVERSE TRANSCRIPTION PROTOCOL

- 7.1** Thaw content of "RT-buffer" and "RT-HAV+HCV+HDV+HGV+HIV primers and dNTP" tubes from Reverse Transcription Reagent Set at room temperature, then vortex thoroughly and spin down drops for 3-5 sec.



The RT buffer can contain the precipitate. Keep it at room temperature with intermittent vortexing until precipitate fully dissolved.

- 7.2** Prepare the mixture of "RT-Buffer", "RT-HAV+HCV+HDV+HGV+HIV primers and dNTPs" and reverse transcriptase (RT-mix ). Add into the one plastic tube:  
 2.0 x (N+1)  $\mu\text{L}$  "RT-Buffer";  
 1.0 x (N+1)  $\mu\text{L}$  "RT-HAV+HCV+HDV+HGV+HIV primers and dNTPs";  
 0.5 x (N+1)  $\mu\text{L}$  reverse transcriptase,  
 N – the number of samples being analyzed, considering "C-".



Reverse transcriptase should be kept out of freezer chamber for as short time as possible.

- 7.3** Vortex RT-mix and spin down drops for 3-5 sec at room temperature.
- 7.4** Add 3.5  $\mu\text{L}$  of the RT-mix to each tube with isolated RNA sample and to "C-" tube. Vortex the tube and spin down the drops.
- 7.5** Place tubes in thermostat and incubate at 40 °C for 30 min, then heat up to 95 °C and leave for 5 min.
- 7.6** Spin the tubes at 16000 x g for 30 sec to collect the drops.

The cDNA preparation is ready for PCR.



The storage of the cDNA preparation is allowed at temperature from minus 18 °C to minus 22 °C for no longer than 1 month.

## 8. PCR PROTOCOL

**8.1** Mark tubes with PCR-mix for each test sample, negative control (C-), positive control (C+) and two tubes for background buffer (applicable to FLASH PCR kits).

For example if you need to test 10 samples, mark 12 tubes (10 for each sample, 1 for “C-“, 1 for “C+“). For FLASH PCR kit mark 14 tubes (10 for each sample, 1 for “C-“, 1 for “C+“ and 2 for background buffer).



Mark only the caps of the tubes when using Rotor-Gene Q Thermal Cycler.

**8.2** Thaw PCR-buffer at the room temperature.

**8.3** Mix the PCR-buffer and TECHNO Taq-polymerase thoroughly (3-5 sec), then spin briefly (1-3 sec) at room temperature.



Hold TECHNO Taq-polimerase at room temperature as short time as possible. The overheating is detrimental to its performance.

**8.4** Prepare the mixture of PCR-buffer and TECHNO Taq-polymerase (Taq-polymerase solution). Add into the one tube:

- 10 x (N+1) µL of PCR-buffer;
- 0.5 x (N+1) µL of TECHNO Taq-polymerase,

N — number of the marked tubes including C-, C+, background tubes.

For example if you need to test 10 samples (12 marked tubes), prepare mixture of PCR-buffer and TECHNO Taq-polymerase for 13 (12+1) tubes: 130 µL PCR-buffer + 6.5 µL TECHNO Taq-polymerase.

**8.5** Vortex the tube with Taq-polymerase solution for 3-5 seconds and spin down drops for 3-5 sec. The maximum storage time for Taq-polymerase solution is 1 hour.

**8.6** Add 10 µL of Taq-polymerase solution into each tube (except background tubes). Add 10 µL of PCR buffer into corresponding tubes (applicable to FLASH PCR kits). Avoid paraffin layer break.

**8.7** Add one drop (~20 µL) of mineral oil into each tube (skip this step when using R3-P609-24/9EU). Close tubes tightly.

**8.8** Vortex the tubes with samples, “C-“ and “C+“ for 3-5 sec and spin down drops for 3-5 sec.

**8.9** Add 5.0 µL of cDNA sample into corresponding tube. Avoid paraffin layer break. Do not add DNA into the “C-“, “C+“ and background (applicable to FLASH PCR kits) tubes. Avoid paraffin layer break.



Open the tube, add DNA sample, then close the tube before proceeding to the next DNA sample to prevent contamination. Use filter tips.

**8.10** Add 5.0 µL of “C-“ which has passed NA isolation stage and reverse transcription reaction into “C-“ and background (applicable to FLASH PCR Detection Kits) tubes. Add 5.0 µL of “C+“ into corresponding tube. Avoid paraffin layer break.

- 8.11** Spin tubes briefly (1-3 sec) at room temperature (not applicable to kits approved for use with Rotor-Gene Q thermal cycler).
- 8.12** Set the tubes to the Thermal Cycler.
- 8.13** Launch the Thermal Cycler software and run PCR according to instructions supplied with device, considering 35 µl reaction mix volume. See tables 4-8 to refer the cycling program and table 9 to refer the detection channels (applicable to Real-Time PCR kits). Using Thercyc cycler you need to choose “Rapid active regulation” regulation algorithm. Amplification programs correspondence with reagent kits is shown in Table 4.
- In case you use DTlite and DTprime Thermal Cyclers, launch the RealTime\_PCR application in “Device operating” mode. Upload ini file before the first run. Add corresponding test in subsequent runs. Specify the number and identifiers of samples. Define position of tubes in software interface according to position they were set in thermal unit. Run RT-PCR.

Table 4. The PCR program for Tercyc Conventional PCR Thermal Cycler

Step	For thermal cyclers with active regulation			Number of cycles
	Temperature, °C	Time		
		min	sec	
1	94	5	0	1
2	94	0	10	50
	58	0	25	
	64	0	15	
3	10	...	...	Holding



When working with FLASH PCR Kit once prepared and amplified background tubes may be used many times at each PCR results detection with reaction tubes from the same lot. Background tubes can be stored at temperature between 2 °C and 8 °C and out of light for 1 month. During the detection procedure background tubes must be room temperature, for that take out them from refrigerator 1 hour before detection.

Table 5. The PCR program for DTlite and DTprime Thermal Cyclers

Step	Temperature, °C	Min.	Sec.	Number of cycles	Optical measurement	Type of the step
1	94	5	00	1		Cycle
2	94	0	10	50		Cycle
	58	0	25		v	
	64	0	15			
3	10 <sup>1</sup>	...	...	Holding		Holding

<sup>1</sup> – holding at 25°C is allowed

Table 6. The PCR program for iCycler iQ5 thermal cyclers (with persistent well factor)

Cycle	Repeats	Step	Dwell time	Setpoint, °C	PCR/Melt Data Acquisition
1	1				
		1	05:00	94,0	
2	50				
		1	00:10	94,0	
		2	00:25	58,0	Real Time
		3	00:15	64,0	
3	1				
		1	01:00	10,0	

Table 7. The PCR program for iCycler iQ thermal cyclers (Bio-Rad Laboratories)

Cycle	Repeats	Step	Dwell time	Setpoint, °C	PCR/Melt Data Acquisition
dynamicwf.tmo program					
1	1				
		1	00:30	80,0	
		2	05:00	94,0	
2	10				
		1	00:20	94,0	
		2	00:20	58,0	
		3	00:10	64,0	
3	2				
		1	00:20	85,0	Real Time
PCR program					
4	40				
		1	00:10	94,0	
		2	00:10	58,0	
		3	00:30	58,0	Real Time
		4	00:20	64,0	
5		...	...	10,0	Storage

Table 8. The PCR program for Rotor-Gene Q Thermal Cyclers

Cycling	Temperature	Hold Time	Cycle Repeats
Cycling	80 °C	300 sec	1 time
Cycling 2	94 °C	10 sec	50 times
	54 °C	5 sec	
	58 °C <sup>3</sup>	25 sec	
	64 °C	20 sec	

<sup>3</sup> - take the measurement

Table 9. Detection channels

	Specific product	IC
DTprime, DTlite and IQ5	Fam	Hex
Rotor-Gene Q	Green	Yellow

## 9. CONTROLS

Table 10.

Control	The controlled step	Result	Interpretation
C+	PCR	Specific signal on Fam/Green is present and the value fits in the range	Valid
		Specific signal on Fam/Green is absent	Invalid*
C-	PCR and RNA extraction	Specific signal on Fam/Green is present	Invalid**
		Specific signal on Fam/Green is absent	Valid
IC	PCR and RNA extraction	Specific signal on Hex/Yellow is present and the Cp value fits in the range	Valid
		Specific signal on Hex/Yellow is absent	Invalid***

The sample is considered positive if the signal for specific cDNA is present. The signal for IC could be absent in samples with high concentration of specific cDNA due to competitive priming.

The sample is considered negative if the signal for specific cDNA is absent and for IC is present.

\* - if the signal for C+ is absent or the value does not fit in the range, all tests of current batch must be repeated.

\*\* - If the signal for C- is present or the Cp value does not fit in the range, the result can not be reported to the patient and the test should be repeated for the sample.

## 10. DATA ANALYSIS

In case of using DNA-Technology made Real-Time PCR Thermal Cyclers or Fluorescence Readers the analysis performed automatically. In all other cases the analysis is based on the presence or absence of specific signal. The controls should be also considered to exclude false positive and false negative results (see p. 7 of the current manual). The cutoff Ct values for Rotor-Gene Q thermal cycler are 40 (specific product), 33 (C+) and 36 (IC). The result characterized by Ct above this value should be considered doubtful and the whole assay should be repeated.

The interpretation should be performed in accordance with table 11.

Table 11

HIV FLASH PCR detection Kit	Test samples		Interpretation
	HIV REAL-TIME PCR detection Kit		
	Fam (Green for Rotor-Gene Q)	Hex (Yellow for Rotor-Gene Q) <sup>4</sup>	
Test samples			
“+”	Cp/Ct is specified, for Rotor-Gene Q Ct<40	Not considered	Positive
“-”	Cp/Ct is not specified (for iQ N/A)	Cp/Ct 29-34, for Rotor-Gene Q Ct<36	Negative
“uncertain”	Cp/Ct is not specified (for iQ N/A)	Cp/Ct is not specified (for iQ5 N/A), for Rotor-Gene Cp/Ct is not specified or Ct>36	Uncertain
<b>C+</b> <sup>5</sup>			
“+”	Cp/Ct<33	Not considered	Positive
<b>C-</b>			
“-”	Cp/Ct is not specified (for iQ N/A)	Cp/Ct 29-34, for Rotor-Gene Q Ct<36	Negative

<sup>4</sup> - if the Cp (Ct) value for HEX channel exceeds the indicated range the result should be considered uncertain.

<sup>5</sup> - if the Cp (Ct) value for controls exceeds the indicated range the result should be considered uncertain.

## 11. TROUBLESHOOTING

Table 12.

	Specific signal +	Specific signal -	Possible cause	Solution
C+	-	-	Operation error PCR inhibition Violation of storage and handling requirements	Repeat whole test  Dispose current batch
C-	+	+	Contamination	Dispose current batch  Perform decontamination procedures
RNA-IC		-	PCR inhibition	Repeat whole test

If you face to any undescribed issues contact our representative

## 12. STORAGE AND HANDLING REQUIREMENTS

Shelf life - 12 months if all the conditions of transportation, storage and operation are met.

All components of the **HIV PCR detection Kit** except PCR-mix and C+ must be stored at minus 20 °C over the storage period. The PCR-buffer and mineral oil can be stored at temperature between 2°C and 8 °C.

The PCR-mix, C+ and *PREP-NA* DNA/RNA Extraction Kit must be stored at temperature between 2°C and 8 °C over the storage period.

Transportation can be held by all types of roofed transport with adherence to above mentioned temperature requirements.

An expired **HIV PCR detection Kit** must not be used.

We strongly recommend following the instructions to get robust and reliable results.

The conformity of the **HIV PCR detection Kit** to the prescribed technical requirements is subject to compliance of storage, carriage and handling conditions recommended by manufacturer.

Contact our customer service by quality issues of the **HIV PCR detection Kit**:

“DNA-Technology” LLC, 117587, Russia, Moscow, int. ter. Municipal District Chertanovo Severnoye, Varshavskoye shosse, 125 Zh, building 5, floor 1, office 12.

Phone/Fax: +7(495)6401771

e-mail: [help@dna-technology.ru](mailto:help@dna-technology.ru), [www.dna-technology.ru](http://www.dna-technology.ru)



### 13. SPECIFICATIONS

- a. **Analytical specificity:** the **HIV PCR detection Kit** allows detection of HIV-1. The samples containing HIV-1 will be defined as positive. The samples not containing HIV-1 will be defined as negative.
- b. **Sensitivity:** not less than 200 copies of HIV DNA per 1 mL of blood plasma.
- c. **Diagnostic sensitivity:** 99,5%.
- d. **Diagnostic specificity:** 100%.



The claimed specifications are guaranteed when DNA extraction is performed with *PREP-NA* DNA/RNA Extraction Kit.

### 14. QUALITY CONTROL

“DNA-Technology, Research&Production” LLC declares that the above mentioned products meet the provision of the Council Directive 98/79/EC for *In Vitro* Diagnostic Medical Devices. The quality control procedures performed in accordance with ISO 9001:2008 and ISO 13485:2003

## 15. KEY TO SYMBOLS



Authorized representative in EU



Caution



Consult instructions for use



Date of manufacture



Expiration date



For research use only



Batch code



Version



Upper limit of temperature



Manufacturer



Negative control



Positive control



Catalogue number



Sufficient for



Temperature limitation

**Nomogram and formula for calculation of relative centrifugal force (RCF) in the speed of rotation (RPM) depending of the rotor diameter**

